

Product Datasheet

Human SIRP gamma / CD172g Protein (orb257847)

Description Signal-regulatory protein gamma (SIRPG or SIRP-gamma) is also known as CD172 antigen-like family mem...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target SIRP gamma

Preservatives 50 mM Tris, 100 mM Glycine, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 63.7 kDa. The protein migrates as 67-85 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH64532


Purity 95%

Source Human SIRP gamma, Fc Tag (orb257847) is expressed from human 293 cells (HEK293). It contains AA Glu 29 - Ser 364 (Accession # AAH64532).

MW 63.7 kDa

NCBI [AAH64532](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human SIRP gamma pr...