


Product Datasheet

Human TIMP1 Protein (orb257869)

Description	TIMP metalloproteinase inhibitor 1 is also known as TIMP1, a tissue inhibitor of metalloproteinases, ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIMP-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 21.5 kDa. The protein migrates as 30-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH00866
Purity	95%
Source	Human TIMP1, His Tag (orb257869) is expressed from human 293 cells (HEK293). It contains AA Cys 24 - Ala 207 (Accession # P01033-1).
MW	21.5 kDa
NCBI	AAH00866
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human TIMP-1 protei...