


Product Datasheet

Human TIMP-2 Protein (orb257870)

Description	TIMP metalloproteinase inhibitor 2 is also known as TIMP2, which belongs to the protease inhibitor I3...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIMP-2
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-8×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 22.6 kDa. The protein migrates as 23-26 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_003246.1
Purity	95%
Source	Human TIMP-2, His Tag (orb257870) is expressed from human 293 cells (HEK293). It contains AA Cys 27 - Pro 220 (Accession # NP_003246.1).
MW	22.6 kDa
NCBI	NP_003246.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TIMP-
2 protei...