


Product Datasheet

Human TIM-3 / HAVCR2 Protein (orb257872)

Description	Human TIM-3, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 200 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIM-3
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 21.8 kDa. The protein migrates as 44-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_116171
Purity	90%
Source	Human TIM-3, His Tag (orb257872) is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 200 (Accession # Q8TDQ0-1).
MW	21.8 kDa
Uniprot ID	Q8TDQ0
NCBI	NP_116171
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human TIM-
3 protein