

Product Datasheet

Human Osteoprotegerin / TNFRSF11B Protein (orb257888)

Description

Human Osteoprotegerin, His Tag is expressed


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Osteoprotegerin

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 44.4 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002537

Purity

95%

Source


Human Osteoprotegerin, His Tag (orb257888) is expressed from human 293 cells (HEK293). It contains AA Glu 22 - Leu 401 (Accession # NP_002537).

MW

44.4 kDa

Uniprot ID
O00300
NCBI
NP_002537
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human OPG protein