

## Product Datasheet

### Human TRAIL R1 / DR4 / TNFRSF10A Protein (orb257894)

**Description**

Human TRAIL R1, His Tag is expressed from


**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

TRAIL R1

**Preservatives**

PBS, pH7.4

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.8 kDa. The protein migrates as 20-23 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

NP\_003835.3

**Purity**

90%

**Source**

Human TRAIL R1, His Tag (orb257894) is expressed from human 293 cells (HEK293). It contains AA Ala 24 - Asn 239 (Accession # O00220-1).

**MW**

24.8 kDa

**Uniprot ID**
**O00220**
**NCBI**
**NP\_003835.3**
**Expiration Date**

6 months from date of receipt.

SDS-PAGE  
analysis of  
Human  
TRAIL R1  
prot...