


Product Datasheet

Human TYRO3 / Dtk Protein (orb257908)

Description	Human TYRO3, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ala 41 - Ser 428 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TYRO3
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 68.2 kDa. The protein migrates as 86-106 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH51756.1 (Q261H)
Purity	95%
Source	Human TYRO3, Fc Tag (orb257908) is expressed from human 293 cells (HEK293). It contains AA Ala 41 - Ser 428 (Accession # AAH51756.1 (Q261H)).
MW	68.2 kDa
Uniprot ID	Q06418
NCBI	AAH51756.1, Q261
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
TYRO3
protein