


Product Datasheet

Human UCH-L1 Protein (orb257912)

Description	Ubiquitin carboxyl-terminal hydrolase isozyme L1 (UCHL1) ,a member of the peptidase C12 family, is a...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	UCH-L1
Preservatives	50 mM Tris, 150 mM NaCl, pH8.0
Form/Appearance	Powder
Storage	-20°C
Tag	N-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 25.7 kDa. The protein migrates as 27-28 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_004172
Purity	95%
Source	Human UCH-L1, His Tag (orb257912) is expressed from E. coli cells. It contains AA Gln 2 - Ala 223 (Accession # NP_004172).
MW	25.7 kDa
NCBI	NP_004172
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
PGP9.5
protei...