


Product Datasheet

Human uPAR / PLAUR Protein (orb257913)

Description	Human uPAR, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 23 - Arg 303 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	uPAR
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-7×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 32.3 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH02788
Purity	92%
Source	Human uPAR, His Tag (orb257913) is expressed from human 293 cells (HEK293). It contains AA Leu 23 - Arg 303 (Accession # AAH02788).
MW	32.3 kDa
Uniprot ID	Q03405
NCBI	AAH02788
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
uPAR
protein