

Product Datasheet

Human VCAM-1 / CD106 Protein (orb257916)

Description Human VCAM-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Phe 25 - Pro 697 (A...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target VCAM-1

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 74.9 kDa. The protein migrates as 85-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_001069

Purity 95%


Source Human VCAM-1, His Tag (orb257916) is expressed from human 293 cells (HEK293). It contains AA Phe 25 - Pro 697 (Accession # NP_001069).

MW 74.9 kDa

Uniprot ID **P19320**

NCBI **NP_001069**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
CD106
protein