



Product Datasheet

Human VEGF-D Protein (orb257928)

Description	Vascular endothelial growth factor D (VEGF-D) is also known as C-fos induced growth factor (FIGF), w...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	VEGF-D
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-9×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 13.4 kDa. The protein migrates as 19-22 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human VEGF-D, His Tag (orb257928) is expressed from human 293 cells (HEK293). It contains AA Phe 93 - Ser 201 (Accession # AAH27948).
MW	13.4 kDa
NCBI	AAH27948
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human VEGF-D protei...



Human VEGF R3, Fc Tag (Cat. No. orb18510...