

Product Datasheet

Human VSIG2 Protein (orb257930)

Description V-set and immunoglobulin domain-containing protein 2 (VSIG2) is also known as Cortical thymocyte-like...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target VSIG2

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-9×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.4 kDa. The protein migrates as 28-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH07313


Purity 95%

Source Human VSIG2, His Tag (orb257930) is expressed from human 293 cells (HEK293). It contains AA Val 24 - Ala 243 (Accession # AAH07313).

MW 24.4 kDa

NCBI [AAH07313](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
VSIG2
protein