

Product Datasheet

Human VSIG4 Protein (orb257931)

Description V-set and immunoglobulin domain-containing protein 4 (VSIG4) is also known as Protein Z39Ig, is a ty...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target VSIG4

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-7×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 30.1 kDa. The protein migrates as 37-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence AAH10525


Purity 95%

Source Human VSIG4, His Tag (orb257931) is expressed from human 293 cells (HEK293). It contains AA Arg 20 - Pro 283 (Accession # AAH10525).

MW 30.1 kDa

NCBI **AAH10525**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
VSIG4
protein