


Product Datasheet

**Ebolavirus EBOV (subtype Zaire, strain Kikwit-95)
Envelope Glycoprotein 1 (GP1) Protein
(orb257933)**

Description

Ebolavirus EBOV (subtype Zaire, strain



SDS-PAGE
analysis of
Viral
Glycoprotein

...

	from human 293 cells (HE...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Glycoprotein / GP (virus)
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-8×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 52.0 kDa. The protein migrates as 95-120 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	AAQ55048.1
Purity	95%
Source	Ebolavirus EBOV (subtype Zaire, strain Kikwit-95) GP1, His Tag (orb257933) is expressed from human 293 cells (HEK293). It contains AA Ile 33 - Arg 501 (Accession # AAQ55048.1).
MW	52.0 kDa
Uniprot ID	P87666
NCBI	AAQ55048.1
Expiration Date	6 months from date of receipt.