


Product Datasheet

Mouse ADAM17 / TACE / CD156b Protein (orb257936)

Description	Mouse ADAM17, His Tag is expressed from human 293 cells (HEK293). It contains AA Pro 18 - Asp 563 (A...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ADAM17
Preservatives	50 mM MES, 150 mM NaCl, pH6.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 63.4 kDa. The protein migrates as 19 kDa and 20 kDa (Propeptide), 55-66 kDa (Mature-ADAM17) and 75-94 kDa (Pro-ADAM17) under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Pro 18 - Asp 563
Purity	90%
Source	Mouse ADAM17, His Tag (orb257936) is expressed from human 293 cells (HEK293). It contains AA Pro 18 - Asp 563 (Accession # Q9Z0F8).
MW	63.4 kDa
Uniprot ID	Q9Z0F8
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse
ADAM17
protei...