

Product Datasheet

Human HABP1 / C1QBP Protein (orb257938)

Description Complement component 1 Q subcomponent-binding protein, mitochondrial (C1QBP), a member of the MAM33 ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target HABP1

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-6×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.6 kDa. The protein migrates as 28-31 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence NP_001203.1


Purity 95%

Source Human HABP1, His Tag (orb257938) is expressed from E.coli cells. It contains AA His 75 - Gln 282 (Accession # Q07021-1).

MW 24.6 kDa

NCBI [NP_001203.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
HABP1
protein