


Product Datasheet

Human DLL1 / Delta1 Protein (orb257943)

Description	Human DLL1, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 18 - Gly 540 (Acc...
Reactivity	Human
Endotoxins	0.01 EU per µg
Conjugation	Unconjugated
Target	DLL1
Preservatives	50 mM Tris, 150 mM NaCl, pH 7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.4 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_005609
Purity	95%
Source	Human DLL1, His Tag (orb257943) is expressed from human 293 cells (HEK293). It contains AA Gln 18 - Gly 540 (Accession # NP_005609).
MW	58.4 kDa
Uniprot ID	O00548
NCBI	NP_005609
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
DLL1
protein