

Product Datasheet

Human DNAM-1 / CD226 Protein (orb257944)

Description Human DNAM-1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 19 - Asn 247 (Ac...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target DNAM-1

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 52.6 kDa. The protein migrates as 66-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_006557

Purity 95%


Source Human DNAM-1, Fc Tag (orb257944) is expressed from human 293 cells (HEK293). It contains AA Glu 19 - Asn 247 (Accession # NP_006557).

MW 52.6 kDa

Uniprot ID [Q15762](#)

NCBI [NP_006557](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
DNAM-1
protei...