


Product Datasheet

Human FKBP12 / FKBP1A Protein (orb257947)

Description	Peptidyl-prolyl cis-trans isomerase FKBP1A, a member of the FKBP-type PPlase family and FKBP1 subfam...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FKBP12
Preservatives	50 mM HEPES, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 12.8 kDa. The protein migrates as 13 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_000792.1
Purity	95%
Source	Human FKBP12, His Tag (orb257947) is expressed from E.coli cells. It contains AA Met 1 - Glu 108 (Accession # P62942-1).
MW	12.8 kDa
NCBI	NP_000792.1
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Human
FKBP12
protei...