

Product Datasheet

Human Dkk-3 (R335G) Protein (orb257976)

Description

Human Dkk-3 (R335G), His Tag (DK3-H5220) is


Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Dkk-3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-6×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.0 kDa. The protein migrates as 50-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_001018067.1 (R335G)

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source


Human Dkk-3 (R335G), His Tag (orb257976) is expressed from human 293 cells (HEK293). It contains AA Ala 22 - Ile 350 (Accession # Q9UBP4-1 (R335G)).

MW

37.0 kDa

Uniprot ID
[Q9UBP4](#)
NCBI
[R335, NP_001018067.1](#)
Expiration Date

6 months from date of receipt.



SDS-PAGE analysis of Human Dkk-3 protein