



Product Datasheet

NSUN7 antibody (orb258175)

Description	Rabbit polyclonal antibody to NSUN7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 19-199aa of human NSUN7 protein
Target	NSUN7
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 lysate using NSUN7 antibody.



Immunohistochemical staining of fetal colon tissue using NSUN7 antibody.