

Product Datasheet

MRPL28 antibody (orb258189)

Description	Rabbit polyclonal antibody to MRPL28
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 91-245aa of human MRPL28 protein
Target	MRPL28
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	30kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.