



Product Datasheet

ST8SIA5 antibody (orb25821)


Description	ST8SIA5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Alpha-2,8-sialyltransferase 8E protein (101-376AA)
Target	ST8SIA5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15466
Dilution Range	WB:1:500-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human prostate tissue.



Immunohistochemical staining of human liver tissue.



Western blot analysis of U87 whole cell lysate.