

Product Datasheet

XPNPEP3 antibody (orb258313)

Description	Rabbit polyclonal antibody to XPNPEP3.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 42-261aa of human XPNPEP3 protein
Target	XPNPEP3
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Concentration	0.5 mg/ml
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Dilution Range	WB: 1:500-1:2,000. IHC: 1:100-1:200. ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.