

Product Datasheet

SNRNP25 antibody (orb258504)

Description	Rabbit polyclonal antibody to SNRNP25.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant fragment corresponding to 1-132aa of human SNRNP25 protein
Target	SNRNP25
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Concentration	0.5 mg/ml
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	15KD
Expiration Date	12 months from date of receipt.