


Product Datasheet

AVL9 antibody (orb258611)

Description	Rabbit polyclonal antibody to AVL9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 2-170aa of human AVL9 protein
Target	AVL9
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles
Note	For research use only
Application notes	Reconstitute in 100µl of sterile H ₂ O with 50% glycerol
Isotype	IgG
Clonality	Polyclonal
Purity	Purified by antigen affinity column
MW	38 kDa
Dilution Range	WB: 1:200-1:1,000, IHC: 1:100 - 1:200, ELISA 1:20,000 - 80,000 - Optimal dilutions/concentrations should be determined by the user
Expiration Date	12 months from date of receipt.



Western blot against recombinant protein...