

Product Datasheet

PPP4R4 antibody (orb258670)

Description	Rabbit polyclonal antibody to PPP4R4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 85-266aa of human PPP4R4 protein
Target	PPP4R4
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	recombinant protein using :45kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.