

Product Datasheet

Tbx15 antibody (orb25883)

Description	Tbx15 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse T-box transcription factor TBX15 protein (1-120AA)
Target	Tbx15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	070306
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.