


Product Datasheet

AURKB antibody (orb259058)

Description	Goat polyclonal antibody to AURKB
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Aurora Kinase B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in Human Thymus lysates and approx. 35kDa in Human Lymph Node, Spleen, Tonsil and Testis lysates (calculated MW of 39.3kDa according to NP_004208.2). Recommended concentration: 1-3µg/ml.
Protein Sequence	YKELQKSCTFDEQ
Clonality	Polyclonal
MW	39.3
Entrez	9212
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Thymus (Lane 1)...