


Product Datasheet

IL17C antibody (orb259059)

Description	Goat polyclonal antibody to IL17C
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL17C (aa160-172)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 20kDa band observed in Human Lymph Node lysates (calculated MW of 21.8kDa according to NP_037410.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	RRPCSRDGSGLPT
Clonality	Polyclonal
MW	21.8
Entrez	27189
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lymph Nod...