


Product Datasheet

IL17D antibody (orb259061)

Description	Goat polyclonal antibody to IL17D
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL-17D (aa55-65)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 20 and 23kDa bands observed in Human Brain (Frontal Cortex) lysates (calculated MW of 21.9kDa according to NP_612141.1). The observed 23kDa band corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.
Protein Sequence	HHTLQLGPREQ
Clonality	Polyclonal
MW	21.9
Entrez	53342
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Frontal C...