

Product Datasheet

TM4SF5 antibody (Biotin) (orb25929)

Description	TM4SF5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Transmembrane 4 L6 family member 5 protein (112-157AA)
Target	TM4SF5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14894
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.