

Product Datasheet

Adiponectin Receptor 1/ADIPOR1 Antibody (orb259595)

Description

Adiponectin Receptor 1/ADIPOR1 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ICC, IF, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the N-terminus of human ADIPOR1 (51-78aa EQTCPVPQEEEEVRVLTLPQAHHAME), different from the related mouse sequence by two amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for ADIPOR1 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

42616 MW

Uniprot ID

[Q96A54](#)


Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Rat
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human


Expiration Date

12 months from date of receipt.

1 2 3 4



WB analysis of ADIPOR1 using anti-ADIPOR1...



IF analysis of ADIPOR1 using anti-ADIPOR1...