

Product Datasheet

APH1a Antibody (orb259598)

Description	APH1a Antibody	116KD 1 2 3 4 5
Species/Host	Rabbit	58KD—
Reactivity	Human, Mouse	40KD—
Conjugation	Unconjugated	29KD—
Tested	WB	20KD—
Applications		14KD—
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human APH1a (236-265aa LRSIQRSLLCRRQEDSRVMVYSALRIPPED), different from the related mouse sequence by one amino acid.	WB analysis of APH1a using anti-APH1a an...
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	WB: The detection limit for APH1a is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	28996 MW	
Uniprot ID	Q96BI3	
Dilution Range	Western blot, 0.1-0.5µg/ml, Human, Mouse	
Expiration Date	12 months from date of receipt.	