

Product Datasheet

SNAI3 Antibody (orb259617)

Description		100KD —
Species/Host	Rabbit	55KD —
Reactivity	Human	35KD — 25KD — 15KD —
Conjugation	Unconjugated	
Tested Applications	WB	WB analysis of SNAI3 using anti-SNAI3 an...
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SNAI3 (155-178aa FHCHKPYHTLAGLARHRQLHCHLQ), different from the related mouse sequence by two amino acids.	
Form/Appearance	Lyophilized	
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Application notes	WB: The detection limit for SNAI3 is approximately 0.1ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	32474 MW	
Uniprot ID	Q3KNW1	
Dilution Range	Western blot, 0.1-0.5µg/ml, Human	
Expiration Date	12 months from date of receipt.	