

Product Datasheet

Osteocalcin antibody (orb259644)

Description

Osteocalcin antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Osteocalcin. Please contact us for the exact immunogen sequence. The peptide is available as orb374879.

Target

Osteocalcin

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

11 kDa

Uniprot ID

P86547, **P02818**, **P04640**

NCBI

954642

Entrez



632

Dilution Range


IHC-P: 1:100-300 (based on 0.5 mg/ml), WB: 1:100-500

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human endometrium.



Immunohistochemical staining of paraffin-embedded tissue.



Immunohistochemical staining of human endometrium.