



Product Datasheet

CD44 antibody (orb259645)


Description	CD44 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD44. Please contact us for the exact immunogen sequence. The peptide is available as orb374880.
Target	CD44
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P15379, P16070, P26051
NCBI	001001389, 000601
Entrez	960
Dilution Range	WB: 1:200-2000, IHC-P: 1:100-800
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells (lan...



IHC-P image of human endometrial cancer ...



Immunohistochemical staining of pig hear...