

Product Datasheet

TMEM95 antibody (Biotin) (orb26208)

Description	TMEM95 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Transmembrane protein 95 protein (17-145AA)
Target	TMEM95
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3KNT9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.