

Product Datasheet

TSPAN8 antibody (Biotin) (orb26296)

Description	TSPAN8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tetraspanin-8 protein (110-205AA)
Target	TSPAN8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19075
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.