

Product Datasheet

UBE2V1 antibody (FITC) (orb26338)

Description	UBE2V1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 variant 1 protein (2-147AA)
Target	UBE2V1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13404
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.