

Product Datasheet

ZNF333 antibody (FITC) (orb26514)

Description	ZNF333 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Zinc finger protein 333 protein (21-182AA)
Target	ZNF333
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96JL9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.