

Product Datasheet

Synaptosomal Protein 25kDa (SNAP-25) (a.a. 183-197) cleaved & uncleaved Antibody (orb26634)

Description	Synaptosomal Protein 25kDa (SNAP-25) (a.a. 183-197) cleaved & uncleaved Antibody
Species/Host	Culture
Conjugation	Unconjugated
Tested Applications	EIA, IHC
Target	Synaptosomal Protein 25kDa (SNAP-25) (a.a. 183-197) cleaved & uncleaved
Preservatives	No Preservative
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.