

Product Datasheet

ZFP69B antibody (HRP) (orb26718)

Description	ZFP69B antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Zinc finger protein 69 homolog B protein (15-310AA)
Target	ZFP69B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UJL9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.