

Product Datasheet

alpha 2 Macroglobulin antibody (orb2676)

Description	alpha 2 Macroglobulin antibody
Species/Host	Rabbit
Reactivity	Gallus, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human alpha 2 Macroglobulin (751-850/1474 aa)
Target	alpha 2 Macroglobulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160 kDa
Uniprot ID	P01023
Expiration Date	12 months from date of receipt.