

Product Datasheet

H2AFV antibody (FITC) (orb26833)

Description	H2AFV antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Histone H2A.V protein (1-128AA)
Target	H2AFV
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q71UI9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.