


Product Datasheet

TNFRSF10B antibody (orb26954)

Description	TNFRSF10B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 10B protein (56-210AA)
Target	TNFRSF10B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14763
Dilution Range	IF:1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela whole cell...