

Product Datasheet

HSPA12B antibody (FITC) (orb27121)

Description	HSPA12B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Heat shock 70 kDa protein 12B protein (201-500AA)
Target	HSPA12B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96MM6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.