

Product Datasheet

BHMT2 antibody (FITC) (orb27141)

Description	BHMT2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human S-methylmethionine--homocysteine S-methyltransferase BHMT2 protein (1-200AA)
Target	BHMT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H2M3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.