

Product Datasheet

alpha 1 Antitrypsin antibody (orb27156)

Description	alpha 1 Antitrypsin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Alpha-1-antitrypsin protein (25-418AA)
Target	SERPINA1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01009
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.