

Product Datasheet

Epstein-Barr Virus (LMP) antibody [Out of stock] (orb27217)

Description	Mouse Monoclonal to Epstein-Barr Virus (LMP)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P
Immunogen	Fusion protein of beta-galactosidase and EBV-late membrane antigen
Preservatives	Antibody solution in stabilizing phosphate buffer pH 7.3. Contains 0.09 % sodium azide**. The volume is sufficient for at least 100 immunohistochemical tests (100 µl working solution / test). Use appropriate antibody diluent e.g. BIOLOGO Art .No. PU002.
Concentration	50 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Dilution Range	(liquid conc.) 1:10 - 1:20
Expiration Date	12 months from date of receipt.